


PRODUCT SELECTION CATALOG









RTD™ Wound Dressing is a proprietary highly absorbent antimicrobial polyurethane foam with integrated:

- Methylene Blue (0.25 mg/g)
- Gentian Violet (0.25 mg/g)
- Silver Ion (7 mg/g)
- Capillary-Suction Therapy Foam



RTD™ WOUND DRESSING



PRODUCT SIZE	PRODUCT CODE / PART NUMBER	HCPCS	NUMBER PER BOX	EXUDATE LIGHT TO HEAVY	EXUDATE MODERATE TO HEAVY	TUNNELS / CAVITIES	STACKING	ANTIMICROBIAL / ANTIFUNGAL	ANALGESIC
									

RTD™ WOUND DRESSING- NON ADHESIVE									
2 x 2 x 1/4 inch	RTD76020214	A6209	10		●	●	●	●	●
2 x 2 x 1/8 inch	RTD76020218	A6209	10	●			●	●	●
4 x 4 x 1/4 inch	RTD72040414	A6209	20		●	●	●	●	●
4 x 4 x 1/4 inch	RTD76040414	A6209	10		●	●	●	●	●
4 x 4 x 1/8 inch	RTD76040418	A6209	10	●			●	●	●
4 x 5 x 1/4 inch	RTD76040514	A6210	10		●	●	●	●	●
4 x 5 x 1/8 inch	RTD76040518	A6210	10	●			●	●	●
8 x 16 x 1/4 inch	RTD75081614	A6211	5		●	●	●	●	●

*When using RTD™ for tunneling please ensure that you follow the Instruction for Use and Clinical Recommendations.

***INDICATIONS FOR USE** The RTD™ Wound Dressing is indicated for the treatment of partial to full thickness wounds with moderate to heavy exudate, including pressure ulcers, venous stasis and arterial ulcers, diabetic foot ulcers, skin tears, graft wounds and donor sites, surgical/post-operative wounds, first and second degree burns, lacerations and abrasions and any other wound inflicted by trauma.



PRODUCT SELECTION CATALOG

WHERE TO ORDER:

FOR ASSISTANCE WITH ORDERING
INFORMATION, PLEASE CONTACT
KENERIC HEALTHCARE

EMAIL: CustomerService@KenericHC.com
PHONE: 855.872.2013
FAX: 855.763.0271
WEBSITE: www.KenericHealthcare.com

PLEASE CONTACT US FOR
INFORMATION ON ADDITIONAL
LOCAL DISTRIBUTOR PARTNERS
THAT ARE NOT ON THIS LIST

RTD™ is available nationwide.

DISTRIBUTORS

McKesson Medical Surgical

866.625.2679 | Primary Care
888.822.8111 | Extended Care

Medline

800.633.5463

Mercy Surgical Dressing Group

800.637.2950
412.446.1929 | New Account Setup

Moore Medical

800.234.1464

Owens & Minor

contact local sales

Independence Medical

800.860.8027

Marathon Medical

877.431.4753

Premier Medical Distribution

888.932.0019

MedSupply

800.889.9081

InHome Medical Solutions

877.867.6882

Gericare Medical Supply

800.552.4415

Gulf South Medical Supply

800.347.2456

PSS

800.552.4415

3RD PARTY BILLERS

Prism Medical Products

888.244.6421

Advanced Tissue

877.811.608

McKesson PCS

800.451.6510 National Rehab
888.907.8775 Sterling Medical

Byram Healthcare

877.902.9726

Strive Medical

888.771.9229

Total Wound Solutions

888.858.9988

